

01 ATAC-seq method overview

Grigorios Georgolopoulos

2024-06-23

Table of contents

1	ATAC-seq	1
2	Analysis workflow	2

1 ATAC-seq

ATAC-seq (Assay for Transposase-Accessible Chromatin with high-throughput sequencing) is a chromatin accessibility method for detecting active (accessible) chromatin regions and therefore DNA regulatory elements genome-wide. It utilizes a hyperactive Tn5 transposase to insert sequencing adapters (tagmentation) into open chromatin regions (Figure 1). The tagged DNA fragments are then purified, PCR amplified and sequenced resulting reads that cluster within these regions of increased accessibility. The more reads in a region, the more accessible the chromatin is.

ATAC-seq quickly gained popularity due to the protocol simplicity as it requires no sonication and phenol-chloroform purification like FAIRE-seq, or enzymatic digestion like MNase-seq and DNase I-seq. In addition, ATACseq experiments typically require smaller starting material (~50,000 cells) in comparison to DNase I where typically millions of cells are required. In contrast, ATAC-seq yields lower resolution chromatin accessibility datasets compared to DNase I-seq and is limited by the Tn5 bias where the transposase can show insertion preference into certain DNA sequences




Figure 1: ATACseq method overview from Grandi et al., 2022

Some popular applications of ATAC-seq: 1. Identify open chromatin regions 2. Perform differential accessibility analysis 3. Identify transcription factor binding sites (footprinting) 4. Infer gene activity based on gene body and upstream accessibility

2 Analysis workflow

- FASTQ QC
 - Sequence quality
 - GC content
 - Sequence duplication rate
 - Read-length distribution
 - K-mer enrichment
 - Adapter contamination

- Adapter and low Q bases trimming

Recommended tools: [FASTQC](#), [fastp](#), [Pear](#), [cutadapt](#), [trimmomatic](#)

- Alignment & QC
 - Mapping (BWA/bowtie2)
 - Filtering mapped reads (MAPQ ≥ 30)
 - Check alignment rate
 - Count % duplicate reads and filter out
 - Obtain fragments
 - Check insert sizes
 - Generate coverage tracks (bigWigs)
 - Check TSS enrichment

Recommended tools: [BWA-MEM](#), [bowtie2](#), [samtools](#), [picard](#), [deepTools](#)

- Peak calling
 - Call peaks (MACS2/3)
 - Annotate peaks
 - Motif enrichment

Recommended tools: [MACS](#), [HOMER](#), [GREAT](#), [MEME](#)

- TF Footprinting
 - De novo footprint detection
 - Motif-centring footprinting

Recommended tools: [TOBIAS](#), [ChromVar](#), [HINT-ATAC](#)